



Ehrenstorfer Quality



Smart Solutions™ v59 PharmaVetResiMix Kit

The easiest and quickest way
to calibrate and screen for 59
Veterinary Residue analytes.

Find out more at:
lgcstandards.com

LGC Quality | ISO 17034 | ISO/IEC 17025 | ISO 9001

59 ANALYTES . 30 MINUTES . MAXIMUM RESULTS .

The Dr. Ehrenstorfer Smart Solutions™ **v59 PharmaVetResiMix Kit** - 59 analytes, 4 mixtures, 1 kit. Quick calibration and mass screening for LC/MS.

Improve efficiency with the largest mix of its kind. Screen 59 analytes in one 30-minute run:

- Increase productivity, complete more runs in a shift
- Reduce calibration runs and save preparation time

Free up resources for other tasks with our complete solution. Precisely formulated to mitigate issues with method development:

- Increase efficiency with higher sample throughput
- Save time compared to preparing solutions in-house

Have confidence in your results. Produced in one of our ISO 17034 accredited manufacturing facilities:

- Reduce the likelihood of downtime and reduce issues
- Avoid error-prone preparations



Part No.	Description	UoM	Format	No. of Analytes	Shelf Life
DRE-K30000030	Smart Solutions™ v59 PharmaVetResiMix Kit 1	4 x 1 mL	KIT	59	36 months
DRE-A30000031ME	PharmaVetResiMix 1 Avermectins-v6 10 µg/mL in Methanol	1 mL	Mixture	6	36 months
DRE-A30000032ME	PharmaVetResiMix 2 Benzimidazoles-v14 10 µg/mL in Methanol	1 mL	Mixture	14	36 months
DRE-A30000033ME	PharmaVetResiMix 3 Quinolones-v18 10 µg/mL in Methanol	1 mL	Mixture	18	48 months
DRE-A30000034ME	PharmaVetResiMix 4 Sulfonamides-v21 10 µg/mL in Methanol	1 mL	Mixture	21	60 months

Complimentary mixes also available to support your analysis:

Part No.	Description	UoM	Format	No. of Analytes	Shelf Life
DRE-A30000035ME	PharmaVetResiMix 5 Tetracyclines-v10 10 µg/mL in Methanol	1 mL	Mixture	10	36 months
DRE-A30000036WL	PharmaVetResiMix 6 Betalactams-v23 0.2-2.5 µg/mL in Water:Acetonitrile	1 mL	Mixture	23	4 months

WHAT IS YOUR SMART SOLUTION™?

Challenge us to meet your analyte needs, no matter the complexity. Put us to the test at [lgcstandards.com/custom](https://www.lgcstandards.com/custom)

