



Ehrenstorfer Quality



## SMART SOLUTIONS™ TERPIMIX KIT

TerpiMix: 63 analytes, one  
kit: designed to boost your  
laboratory's productivity.

Request a custom organic solution quote at:  
[lgcstandards.com/custom](https://lgcstandards.com/custom)

LGC Quality | ISO 17034 | ISO/IEC 17025 | ISO 9001



## 63 ANALYTES. 1 TEST. ENSURING HIGH-QUALITY RESULTS.

Terpenes are naturally occurring compounds that contribute to the aromas and flavour of cannabis and food products. They also have therapeutic properties and can support the efficacy of certain cannabinoids. The Dr. Ehrenstorfer Smart Solutions™ TerpiMix Kit is designed for optimum stability and elution in order to support the cannabis and food testing industry.

We are excited to be introducing the NEW Dr. Ehrenstorfer Smart Solutions™ TerpiMix Kit 2500 µg/mL in Isopropanol. A set which includes 63 analytes across three mixtures to enable fast multi-analyte screening in cannabis products and the food testing market. These kits are designed to enable optimal elution to ensure accurate analytical interpretations.

- Achieve precise and accurate results with a mixture specifically formulated to remove unnecessary enantiomer isomers.
- Save money by purchasing as a kit or individually to meet your testing requirement needs. Long shelf life fights against waste allowing you to store unused mixtures for over a year.
- Reduce dilution steps with precisely formulated mixtures.
- 63 analytes allow for easier to understand results using spiked method validation.
- Guaranteed quality and easy to reproduce results due to the TerpiMix's ISO 17034 accreditation.
- Test up to 63 terpene analytes at once using the largest and most stable mix of its kind on the market today.



Part No.	Description	UoM	Format	No. of Analytes	Shelf life
DRE-K500006551P <b>NEW</b>	Smart Solutions™ TerpiMix Kit 2500 µg/mL in Isopropanol	3x1ml	Kit	63	550 days
DRE-A500006521P <b>NEW</b>	TerpiMix 1 2500 µg/mL in Isopropanol	1ml	Mixture	22	550 days
DRE-A500006531P <b>NEW</b>	TerpiMix 2 2500 µg/mL in Isopropanol	1ml	Mixture	22	1,825 days
DRE-A500006541P <b>NEW</b>	TerpiMix 3 2500 µg/mL in Isopropanol	1ml	Mixture	19	1,825 days

### WHAT IS YOUR SMART SOLUTION™?

Challenge us to meet your analyte needs, no matter the complexity. Put us to the test at [lgcstandards.com/custom](https://lgcstandards.com/custom)

