


ColdBlock™ Sample Digestion Technology

Innovative
Very rapid
Accurate
Cost-effective

ColdBlock™ Digestion Technology uses short-wave infrared radiation with an integrated cooling block, setting a new performance standard for sample digestion.




ColdBlock™ Digestion offers:


- Rapid digestions (in minutes)
- Excellent accuracy and reproducibility
- Superior safety – reduces need for hazardous reagents
- Greener process
- Easy to scale up and automate

ColdBlock™ Digestion offers unparalleled economic and operational efficiencies via:

- Broad Applicability; mining, metals and metal alloys, environmental, food, agriculture
- High Versatility
- Accelerated sample throughput
- Operational cost savings
- Simple operation



CB12L – 12-channel unit



CB15S – 15-channel unit



Industrial
VHG | ARMI | MBH
Paragon Scientific

lgcstandards.com/industrial - industrial@lgcgroup.com